
Ch-MBE GUI 与 RHEED 后处理标注

Ch-MBE GUI and RHEED Post-processing Labeling

核心边界 / Core boundary
后处理标注器只读取归档数据，不连接或控制相机、温度计、电源及其他仪器。
The post-processing labeler reads archived data only. It does not connect to or control cameras, pyrometers, power supplies, or other instruments.

快速入口 / Quick entry	双击 / Double-click
Ch-MBE 实时 GUI / Live GUI	Start Ch-MBE Growth Monitor.cmd
RHEED 离线标注器 / Offline labeler	Start RHEED Post-processing Labeler.cmd
安装或更新桌面快捷方式 / Install or refresh shortcuts	Install AI4MBE Desktop Shortcuts.cmd

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Understand the two independent paths

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The live GUI performs acquisition and display. The offline labeler works only on completed archived sessions. Never interpret offline output as a real-time control instruction.

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Reconstruction labels, acquisition-quality QC, and FeSe film quality are three different concepts.

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Launcher logs are under %LOCALAPPDATA%\AI4MBE\LauncherLogs and show the resolved Python and repository paths. Verify them against that workstation's configuration; do not assume a fixed drive letter or add global environment variables ad hoc.

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RHEED Post-processing and Temporal Labeling

离线构建时间轴并标记重构分段 / Build an offline timeline and label temporal segments

模型输出可见：仅为 model-assisted review，不是 blind-gold。
Model outputs are visible: exports are not eligible as blind-gold labels.

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The archive must contain one session_metadata.json, one heartbeat_log.csv, and the frames/ images referenced by heartbeat rows. Do not manually rename extracted frames.

```
growth_session.zip
  session_metadata.json
  heartbeat_log.csv
  frames/
    heartbeat_000001_...bmp
    heartbeat_000002_...bmp
```

预测表与模型说明必须成对 / Predictions and model specs are positional pairs

文件 / File	必须包含 / Must contain	检查 / Check
Prediction CSV	frame_index, heartbeat_idx, elapsed_s, captured_at_utc, capture_sequence, frame_name, frame_sha256	行顺序与保存帧完全一致
Model-spec JSON	key, title, classes, probability_columns, status, provenance	类别顺序与概率列一一对应
Session ZIP	metadata, heartbeat, frames	每个 frame_path 在 ZIP 内存在

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Prediction frame_index may begin at 0 or 1, then must remain contiguous. heartbeat_idx, capture_sequence, and elapsed_s may have gaps. The tool follows actual saved frames and recorded times; it never assumes exact 1 Hz sampling.

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Any mismatch in row count, order, time, sequence, filename, or SHA-256 stops the build. Do not edit inputs to bypass the error.

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The report directory contains interactive_report.html, images/, vendor/, and run_manifest.json. Copy the whole directory when moving or handing it off, and preserve the original prediction CSV and model-spec JSON separately outside the report directory.

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
The desktop UI always refuses a non-empty output directory. Advanced CLI --overwrite may replace only a recognized tool output containing run_manifest.json and interactive_report.html, and it validates the new inputs before replacement. It still refuses filesystem roots, repository/source ancestors, the current directory, and directories containing inputs.

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Label temporal segments like an editing timeline

RHEED reconstruction timeline / RHEED 重构时间轴

Model-assisted review - not blind-gold / 模型辅助审阅 - 非盲标金标准

1

拖动 Frame playhead 或点击任意模型曲线，找到分段的第一张保存帧。点击 Mark In。
Move the Frame playhead or click a model plot to find the first saved frame. Click Mark In.

2

移动到该分段的最后一张保存帧，点击 Mark Out。Out 是包含端点。
Move to the final saved frame in the segment and click Mark Out. Out is inclusive.

3

选择重构标签，填写 Labeler，按需添加 Notes，然后点击 Add segment。
Choose a reconstruction label, enter Labeler, add Notes if useful, then click Add segment.

4

点击已保存分段可修改；Update 保留 annotation_id。相邻分段允许，重叠分段会被拒绝。
Select a saved segment to edit it. Update preserves annotation_id. Adjacent segments are allowed; overlaps are rejected.

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Visible and exported saved-frame ordinals are 1-based. Every endpoint also records heartbeat, capture sequence, UTC, and frame SHA-256, so it remains traceable despite sampling gaps.

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Review, display controls, and drafts

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草稿与编辑 / Draft and editing

- 浏览器按 dataset_id 在 localStorage 中保存草稿；它不是长期备份。
The browser stores a dataset-scoped draft in localStorage; this is not a durable backup.
- 频繁 Export JSON，特别是在长会话、换电脑、换浏览器或清理缓存前。
Export JSON frequently, especially for long sessions and before changing computers, browsers, or cache settings.
- Delete selected 需要确认；Undo 只恢复最近一次标注修改。
Delete selected requires confirmation. Undo restores only the most recent annotation mutation.
- Import JSON 会先严格验证；失败时当前分段集合保持不变。
Import JSON validates before replacement. On failure, the current segment set remains unchanged.

标签解释 / Interpretation

标签描述审阅者在该时间区间看到的主要表面重构类别。unknown 表示无法可靠判断。不要用模型 argmax 自动覆盖人的选择。
Labels describe the dominant surface-reconstruction category seen by the reviewer over that interval. Use unknown when a reliable decision is not possible. Never overwrite the human choice with model argmax.

显示增强不是数据增强 / Display adjustment is not data modification
亮度、对比度和放大只帮助观察。导出会记录显示设置，但不会改写归档图像。
Brightness, contrast, and zoom aid inspection only. Exports record display settings but do not rewrite archived images.

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Export, import, and fail-closed validation

JSON 是 canonical 文件 / JSON is the canonical artifact

Export JSON 用于继续标注、验证和下游处理；Export CSV 便于会议和表格分析，但不能替代 JSON 的完整嵌套溯源。
Use Export JSON for resume, validation, and downstream processing. Export CSV is convenient for meetings and tabular analysis, but it does not replace the complete nested provenance in JSON.

绑定字段 / Binding	作用 / Purpose
dataset_id + source_archive_sha256	绑定原始会话 / Bind the source session
ordered_frame_fingerprint	绑定帧顺序与逐帧 SHA / Bind order and per-frame SHA
model_review_context_fingerprint	绑定审阅时可见的 predictions/specs / Bind visible model context
endpoint provenance	绑定 frame ordinal、heartbeat、UTC、sequence 和 SHA / Bind each boundary
model_outputs_visible=true	披露模型可见 / Disclose model visibility
eligible_for_gold=false	禁止误当 blind-gold / Prevent blind-gold misuse

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验证失败不要绕过 / Never bypass a failure
wrong run、模型上下文变化、端点被改、重叠或 gold 声明都会失败。回到正确的报告和原始导出重新核对。
A wrong run, changed model context, altered endpoint, overlap, or gold claim all fail validation. Return to the correct report and original export and investigate.

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Data safety and scientific boundaries

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Keep the source session ZIP read-only or preserve an immutable copy, and record its SHA-256.
- 2

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Create a dedicated working directory outside the repository. Keep reports, exports, and screenshots there.
- 3

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For handoff, copy the whole report directory and canonical JSON, and separately preserve the original prediction CSV, model-spec JSON, and SHA manifest; never send only the HTML.
- 4

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Run statistics or training-data preparation only after validation, and preserve the original export unchanged.
- 不是 blind-gold / Not blind-gold

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Model plots and lossy review images are visible during labeling, so these labels are model-assisted review only. Formal blind-gold data requires a separate workflow that hides model outputs.
- 不要扩大结论 / Do not over-interpret

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The current classification target is bare STO surface reconstruction. Do not directly interpret it as FeSe film quality, and do not equate reconstruction unknown with QC reject.

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First preserve the exact error text, branch/commit, and selected paths. Do not begin by deleting environments, changing global variables, or force-resetting Git. If AI4MBE live-driver warning appears, record its missing-module list and launcher log. Dummy and unaffected modes remain available, but do not ARM a production session that needs a missing driver until the environment is repaired. This is expected fail-visible behavior, not a prompt for ad-hoc installation.

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Version verification and one-page checklist

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Use the resolved Python and repository paths printed by the launcher log. Read AI4MBE_GUI_PYTHON only when it is already configured on that workstation; do not create it ad hoc. If branch/commit differs from the team-specified version, or git status shows unknown changes, stop and verify.

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Stop and contact the maintainer if instrument state, stale data, version identity, provenance, or validation cannot be explained.

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All interface and RHEED images in this manual are synthetic illustrations and contain no experimental data.